

ADJUSTABLE GUNNEL MOUNT FILET TABLE

P01-2120W & P01-2132W

Installation Instructions

Congratulations on your purchase of a TACO Marine Adjustable Gunnel Mount Filet Table!

The Adjustable Gunnel Mount Filet Table is available in 20" and 32" widths. These tables feature a convenient knife and plier holder, and the adjustable tilt allows for desired position on rod holders of 0°, 15° and 30°. Designed openings around the table surface allow for drainage overboard when cleaning fish. The rugged, easy-grip handle adjusts with a single turn and folds flat for convenient storage. To withstand the harsh marine environment, these tables are made of durable, FDA-approved, UV-stabilized, high-density polyethylene.

Box Includes

- TABLE TOP ASSEMBLY
- ADJUSTABLE MOUNT ASSEMBLY
- 1.25" FLAT, PHILLIPS HEAD BOLTS
- WASHERS
- NUTS

Required*
*not included


- ADJUSTABLE WRENCH
- PHILLPS HEAD SCREW DRIVER

ASSEMBLY

- 1) The tilt mount comes as an entire assembly. It is attached with four (4) flat, Phillips head screws.
- 2) Pass the supplied bolts through the top of the table base. (Make sure the adjustable knob faces the front of the table.)
- 3) Install the supplied washers.
- 4) Install the supplied nuts and finger-tighten them as much as possible.
- 5) Use the adjustable wrench to hold the nuts secure.
- 6) Tighten the screws from the top surface of the table using the Phillips head screwdriver.
 - (Be careful not to over-tighten and strip the bolt threads.)

MAINTENANCE

- Hosing the filet table with fresh water is recommended, especially when used in a saltwater environment.
- For stubborn stains, use the same marine-safe cleanser for boats or fishing gear.
- When storing the table, loosen the adjustable mount assembly.
- Lay it flat on the table with the bracket assembly loose to eliminate unneeded stress.



Explore our entire line of Poly Accessories by visiting our Website **tacomarine.com**.

Also, discover a library of exciting how-to videos and more on our popular YouTube channel **youtube.com/tacomarine**.